

# CELLO HOUSE

Disposable shelters for rats, guinea pigs...

**Disposable enrichment shelters stand as the ultimate solution for the welfare of the rodents and animal technicians.**

**Cost effective**, Cello House does not require any time-consuming and costly cleaning after use

**Clean**, Cello House is made of virgin cellulose fibers

**Certified**, Cello House is screened for contaminants

**Very Light**, handling and change are easy and painless

**Clean and Light Material**, Cello House allows shredding, nibbling, moving the house.

**Smooth surface**, a stack of 50 Cello House is half the height of other molded cellulose shelters.

**Environment friendly**, Cello House is made from renewable material and comply with PEFC standards

**Autoclavable**



SERLAB UK Ltd

Lincoln House - ccfgb

300 High Holborn

London WC1V 7JH

Phone : +44 (0) 207 0923436

Email : sales@serlab.co.uk

[www.serlab.co.uk](http://www.serlab.co.uk)

## DESCRIPTION

Arrow shaped shelter made from densified cellulose :  
185 x185x250mm / Height 90mm

Two holes : one on the top 50mm diam and one on  
the side width 75mm - H 50mm

Net weight : 25gr approx

## PACKAGING

108 pieces per carton

Grossweight per carton : 3kg

48 cartons on a 100x120x195 pallet - Gross weight  
165 kg