

LAB COB 8

DESCRIPTION

Round grains obtained through the mechanical grinding and milling of non chemically treated corn cob.



PHYSICAL PROPERTIES

Origin : France.
Raw material : non chemically treated corn cob
Density : 0.350
Moisture content : 12% max
Particule size : 95% between 2 and 3.15 mm
Absorption 105%

SERLAB UK Ltd

Lincoln House - ccfgb

300 High Holborn

London WC1V 7JH

Phone : +44 (0) 207 0923436

Email : sales@serlab.co.uk

www.serlab.co.uk

PACKAGING

12.5 kg autoclavable polypropylene sack
Every bag labelled with a batch number
55 sacks per shrink wrapped pallet

GLP Certificate
Irradiation : optional